

Product Datasheet

Rat Bcl2 ELISA Kit (orb436867)

Description	Rat B-Cell Leukemia/Lymphoma 2 (Bcl2) Sandwich ELISA Kit
Reactivity	Rat
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P49950
Sensitivity	30.98pg/mL
Entrez	24224
Expiration Date	Please enquire.