

Product Datasheet

Mouse Bcl9 ELISA Kit (orb436868)

Description	Mouse B-Cell CLL/Lymphoma 9 (Bcl9) Sandwich ELISA Kit
Reactivity	Mouse
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q9D219
Sensitivity	27pg/mL
Entrez	77578
Expiration Date	Please enquire.