

Product Datasheet

Mouse BLC1 ELISA Kit (orb436878)

Description	Mouse B-Lymphocyte Chemoattractant 1 (BLC1) Sandwich ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	O55038
Sensitivity	5.9pg/mL
Entrez	55985
Expiration Date	Please enquire.