

Product Datasheet

Human BMP2 ELISA Kit (orb436883)

Description	Human Bone Morphogenetic Protein 2 (BMP2) Sandwich ELISA Kit
Reactivity	Human
Range	46.875-3000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P12643
Sensitivity	16.12pg/mL
Entrez	650
Expiration Date	Please enquire.