

Product Datasheet

Rat BMP2 ELISA Kit (orb436886)

Description	Rat Bone Morphogenetic Protein 2 (BMP2) Sandwich ELISA Kit
Reactivity	Rat
Range	62.5-4000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P49001
Sensitivity	25.7pg/mL
Entrez	29373
Expiration Date	Please enquire.