

Product Datasheet

Canine BMP4 ELISA Kit (orb436889)

Description	Canine Bone Morphogenetic Protein 4 (BMP4) Sandwich ELISA Kit
Reactivity	Canine
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	F1PDL1
Sensitivity	13.7pg/mL
Entrez	490695
Expiration Date	Please enquire.