

Product Datasheet

Chicken BMP4 ELISA Kit (orb436890)

Description	Chicken Bone Morphogenetic Protein 4 (BMP4) Sandwich ELISA Kit
Reactivity	Gallus
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	E1C0T2
Sensitivity	12.6pg/mL
Entrez	396165
Expiration Date	Please enquire.