

Product Datasheet

Mouse BMPR1A ELISA Kit (orb436898)

Description	Mouse Bone Morphogenetic Protein Receptor 1A (BMPR1A) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P36895
Sensitivity	0.061ng/mL
Entrez	12166
Expiration Date	Please enquire.