

Product Datasheet

Rat BMPR1A ELISA Kit (orb436899)

Description	Rat Bone Morphogenetic Protein Receptor 1A (BMPR1A) Sandwich ELISA Kit
Reactivity	Rat
Range	0.781-50ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q78EA7
Sensitivity	0.24ng/mL
Entrez	81507
Expiration Date	Please enquire.