

Product Datasheet

Human C1QTNF1 ELISA Kit (orb436900)

Description	Human C1q And Tumor Necrosis Factor Related Protein 1 (C1QTNF1) Sandwich ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q9BXJ1
Sensitivity	0.57ng/mL
Entrez	114897
Expiration Date	Please enquire.