

Product Datasheet

Bovine CTGF ELISA Kit (orb436905)

Description	Bovine Connective Tissue Growth Factor (CTGF) Sandwich ELISA Kit
Reactivity	Bovine
Range	62.5-4000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	018739
Sensitivity	30.23pg/mL
Entrez	281103
Expiration Date	Please enquire.