

Product Datasheet

CD11c antibody (PerCP) (orb43692)

Description

Mouse monoclonal antibody against CD11c conjugated

Reactivity

Human, Monkey

Conjugation

PerCP

Tested

FC

Applications

Immunogen

Dendritic cells of synovial fluid

Target

CD11c

Preservatives

Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide

Storage

Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.

Note

For research use only

Application notes

Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 10⁶ cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.

Isotype

Mouse IgG1

Clonality

Monoclonal

Purity

Purified antibody is conjugated with activated Peridinin-Chlorophyll Protein (PerCP) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.

Uniprot ID

[P20702](#)

Entrez


[3687](#)

Dilution Range

Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 10⁶ cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.

Expiration Date

12 months from date of receipt.



Flow
cytometry
analysis
(surface
stainin...