

Product Datasheet

Rat FGF10 ELISA Kit (orb436937)

Description	Rat Fibroblast Growth Factor 10 (FGF10) Sandwich ELISA Kit
Reactivity	Rat
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P70492
Sensitivity	33pg/mL
Entrez	25443
Expiration Date	Please enquire.