

Product Datasheet

Rat FGF11 ELISA Kit (orb436940)

Description	Rat Fibroblast Growth Factor 11 (FGF11) Sandwich ELISA Kit
Reactivity	Rat
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q8R5L8
Sensitivity	6.2pg/mL
Entrez	170632
Expiration Date	Please enquire.