

Product Datasheet

Mouse FGF21 ELISA Kit (orb436964)

Description	Mouse Fibroblast Growth Factor 21 (FGF21) Sandwich ELISA Kit
Reactivity	Mouse
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q9JJN1
Sensitivity	39.54pg/mL
Entrez	56636
Expiration Date	Please enquire.