

Product Datasheet

Porcine FGF21 ELISA Kit (orb436965)

Description	Porcine Fibroblast Growth Factor 21 (FGF21) Sandwich ELISA Kit
Reactivity	Porcine
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	C6JWG4
Sensitivity	13.3pg/mL
Entrez	100302504
Expiration Date	Please enquire.