

Product Datasheet

Chicken FGF2 ELISA Kit (orb436971)

Description	Chicken Fibroblast Growth Factor 2, Basic (FGF2) Sandwich ELISA Kit
Reactivity	Gallus
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P48800
Sensitivity	6.3pg/mL
Entrez	396413
Expiration Date	Please enquire.