

Product Datasheet

Bovine FGF9 ELISA Kit (orb436996)

Description	Bovine Fibroblast Growth Factor 9 (FGF9) Sandwich ELISA Kit
Reactivity	Bovine
Range	62.5-4000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	D8KXU9
Sensitivity	27.5pg/mL
Entrez	613731
Expiration Date	Please enquire.