

Product Datasheet

Rat FGF9 ELISA Kit (orb437000)

Description	Rat Fibroblast Growth Factor 9 (FGF9) Sandwich ELISA Kit
Reactivity	Rat
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P36364
Sensitivity	37pg/mL
Entrez	25444
Expiration Date	Please enquire.