

Product Datasheet

Human FGFR2 ELISA Kit (orb437004)

Description	Human Fibroblast Growth Factor Receptor 2 (FGFR2) Sandwich ELISA Kit
Reactivity	Human
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P21802
Sensitivity	12.2pg/mL
Entrez	2263
Expiration Date	Please enquire.