

Product Datasheet

Human FGFR1 ELISA Kit (orb437008)

Description	Human Fibroblast Growth Factor Receptor Like Protein 1 (FGFR1) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q8N441
Sensitivity	0.051ng/mL
Entrez	53834
Expiration Date	Please enquire.