

Product Datasheet

Mouse Glucagon ELISA Kit (orb437010)

Description	Mouse Glucagon (GC) competitive inhibition ELISA Kit
Reactivity	Mouse
Range	23.437-1500pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P55095
Sensitivity	17.34pg/mL
Entrez	14526
Expiration Date	Please enquire.