

Product Datasheet

Rat GDF2 ELISA Kit (orb437015)

Description	Rat Growth Differentiation Factor 2 (GDF2) Sandwich ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	M0RB87
Sensitivity	0.126ng/mL
Entrez	290921
Expiration Date	Please enquire.