

Product Datasheet

Mouse GFER ELISA Kit (orb437016)

Description	Mouse Growth Factor, Augmenter Of Liver Regeneration (GFER) Sandwich ELISA Kit
Reactivity	Mouse
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P56213
Sensitivity	29pg/mL
Entrez	11692
Expiration Date	Please enquire.