

Product Datasheet

Human GITR ELISA Kit (orb437017)

Description	Human Glucocorticoid Induced Tumor Necrosis Factor Receptor (GITR) Sandwich ELISA Kit
Reactivity	Human
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9Y5U5
Sensitivity	5.9pg/mL
Entrez	8784
Expiration Date	Please enquire.