

Product Datasheet

Rat GMCSF ELISA Kit (orb437022)

Description	Rat Colony Stimulating Factor 2, Granulocyte Macrophage (GMCSF) Sandwich ELISA Kit
Reactivity	Rat
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P48750
Sensitivity	6.5pg/mL
Entrez	116630
Expiration Date	Please enquire.