

Product Datasheet

General Hyaluronic Acid ELISA Kit (orb437025)

Description	General Hyaluronic Acid (HA) competitive inhibition ELISA Kit
Reactivity	All
Range	6.25-400ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Sensitivity	2.17ng/mL
Expiration Date	Please enquire.