

Product Datasheet

Human HBEGF ELISA Kit (orb437026)

Description	Human Heparin Binding Epidermal Growth Factor Like Growth Factor (HBEGF) Sandwich ELISA Kit
Reactivity	Human
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q99075
Sensitivity	6.3pg/mL
Entrez	1839
Expiration Date	Please enquire.