

Product Datasheet

Canine HGF ELISA Kit (orb437033)

Description	Canine Hepatocyte Growth Factor (HGF) Sandwich ELISA Kit
Reactivity	Canine
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q867B7
Sensitivity	5.6pg/mL
Entrez	403441
Expiration Date	Please enquire.