


Product Datasheet

Human IFNa ELISA Kit (orb437051)

Description	Human Interferon Alpha (IFNa) Sandwich ELISA Kit
Reactivity	Human
Range	7.8125-500pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P01562
Sensitivity	3.0pg/mL
Entrez	3439
Expiration Date	Please enquire.



Human
Interferon
Alpha (IFNa)
ELISA Kit