

Product Datasheet

Human IFN γ ELISA Kit (orb437072)

Description	Human Interferon Gamma (IFN γ) Sandwich ELISA Kit
Reactivity	Human
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P01579
Sensitivity	11.74pg/mL
Entrez	3458
Expiration Date	Please enquire.