

Product Datasheet

Chicken IGF1 ELISA Kit (orb437081)

Description	Chicken Insulin Like Growth Factor 1 (IGF1) Sandwich ELISA Kit
Reactivity	Gallus
Range	18.75-1200pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P18254
Sensitivity	7.4pg/mL
Entrez	418090
Expiration Date	Please enquire.