

Product Datasheet

Guinea pig IGF1 ELISA Kit (orb437082)

Description	Guinea pig Insulin Like Growth Factor 1 (IGF1) Sandwich ELISA Kit
Reactivity	Guinea pig
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P17647
Sensitivity	0.125ng/mL
Entrez	100135551
Expiration Date	Please enquire.