

Product Datasheet

Mouse IGF1 ELISA Kit (orb437084)

Description	Mouse Insulin Like Growth Factor 1 (IGF1) Sandwich ELISA Kit
Reactivity	Mouse
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P05017
Sensitivity	29pg/mL
Entrez	16000
Expiration Date	Please enquire.