

Product Datasheet

Rat IGF1R ELISA Kit (orb437089)

Description	Rat Insulin Like Growth Factor 1 Receptor (IGF1R) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P24062
Sensitivity	0.062ng/mL
Entrez	25718
Expiration Date	Please enquire.