

Product Datasheet

Bovine IGF2 ELISA Kit (orb437090)

Description	Bovine Insulin Like Growth Factor 2 (IGF2) Sandwich ELISA Kit
Reactivity	Bovine
Range	3.125-200ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P07456
Sensitivity	1.44ng/mL
Entrez	281240
Expiration Date	Please enquire.