

Product Datasheet

Rabbit IGF2 ELISA Kit (orb437096)

Description	Rabbit Insulin Like Growth Factor 2 (IGF2) Sandwich ELISA Kit
Reactivity	Rabbit
Range	4.687-300ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	B7NZU3
Sensitivity	1.5ng/mL
Entrez	100328792
Expiration Date	Please enquire.