

Product Datasheet

Human IGFALS ELISA Kit (orb437098)

Description	Human Insulin Like Growth Factor Binding Protein, Acid Labile Subunit (IGFALS) Sandwich ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P35858
Sensitivity	0.65ng/mL
Entrez	3483
Expiration Date	Please enquire.