

Product Datasheet

Porcine IGFBP1 ELISA Kit (orb437103)

Description	Porcine Insulin Like Growth Factor Binding Protein 1 (IGFBP1) Sandwich ELISA Kit
Reactivity	Porcine
Range	3.125-200ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	E0A8P1
Sensitivity	1.32ng/mL
Entrez	397270
Expiration Date	Please enquire.