

Product Datasheet

Bovine IGFBP2 ELISA Kit (orb437105)

Description	Bovine Insulin Like Growth Factor Binding Protein 2 (IGFBP2) Sandwich ELISA Kit
Reactivity	Bovine
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P13384
Sensitivity	33pg/mL
Entrez	282260
Expiration Date	Please enquire.