

Product Datasheet

Rat IGFBP2 ELISA Kit (orb437110)

Description	Rat Insulin Like Growth Factor Binding Protein 2 (IGFBP2) Sandwich ELISA Kit
Reactivity	Rat
Range	3.125-200ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P12843
Sensitivity	1.34ng/mL
Entrez	25662
Expiration Date	Please enquire.