

Product Datasheet

Bovine IGFBP3 ELISA Kit (orb437111)

Description	Bovine Insulin Like Growth Factor Binding Protein 3 (IGFBP3) Sandwich ELISA Kit
Reactivity	Bovine
Range	62.5-4000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P20959
Sensitivity	23.5pg/mL
Entrez	282261
Expiration Date	Please enquire.