

Product Datasheet

Bovine IGFBP4 ELISA Kit (orb437116)

Description	Bovine Insulin Like Growth Factor Binding Protein 4 (IGFBP4) Sandwich ELISA Kit
Reactivity	Bovine
Range	93.75-6000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q05716
Sensitivity	45.87pg/mL
Entrez	282262
Expiration Date	Please enquire.