

Product Datasheet

Human IGFBP4 ELISA Kit (orb437117)

Description	Human Insulin Like Growth Factor Binding Protein 4 (IGFBP4) Sandwich ELISA Kit
Reactivity	Human
Range	3.125-200ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P22692
Sensitivity	1.13ng/mL
Entrez	3487
Expiration Date	Please enquire.