

Product Datasheet

Porcine IGFBP4 ELISA Kit (orb437119)

Description	Porcine Insulin Like Growth Factor Binding Protein 4 (IGFBP4) Sandwich ELISA Kit
Reactivity	Porcine
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P24854
Sensitivity	0.098ng/mL
Entrez	100144490
Expiration Date	Please enquire.