



Product Datasheet

Cytokeratin 8 antibody (orb43712)


Description	Mouse Monoclonal to Cytokeratin 8.
Reactivity	Bovine, Human, Porcine, Rabbit, Sheep
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-P, IP, WB
Immunogen	Cytoskeleton preparation from HeLa human cervix carcinoma cell line.
Target	Cytokeratin 8
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Immunohistochemistry (paraffin sections): Recommended dilution: 10 µg/ml; positive tissue: prostate. Flow cytometry: Recommended dilution: 1-5 µg/ml, intracellular staining.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	Q7L4M3
Entrez	3856
Dilution Range	Immunohistochemistry (paraffin sections): Recommended dilution: 10 µg/ml; positive tissue: prostate. Flow cytometry: Recommended dilution: 1-5 µg/ml, intracellular staining.
Expiration Date	12 months from date of receipt.



Western blotting detection of cytokerati...



Immunocytochemistry staining of cytokerati...



Immunohistochemistry staining of human p...