

Product Datasheet

Rat IGFBP6 ELISA Kit (orb437127)

Description	Rat Insulin Like Growth Factor Binding Protein 6 (IGFBP6) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P35572
Sensitivity	0.064ng/mL
Entrez	25641
Expiration Date	Please enquire.