

Product Datasheet

Mouse IL10 ELISA Kit (orb437137)

Description	Mouse Interleukin 10 (IL10) Sandwich ELISA Kit
Reactivity	Mouse
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P18893
Sensitivity	34.5pg/mL
Entrez	16153
Expiration Date	Please enquire.