

Product Datasheet

Rat IL10 ELISA Kit (orb437139)

Description	Rat Interleukin 10 (IL10) Sandwich ELISA Kit
Reactivity	Rat
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P29456
Sensitivity	27pg/mL
Entrez	25325
Expiration Date	Please enquire.