

Product Datasheet

Mouse IL10Ra ELISA Kit (orb437141)

Description	Mouse Interleukin 10 Receptor Alpha (IL10Ra) Sandwich ELISA Kit
Reactivity	Mouse
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q61727
Sensitivity	29pg/mL
Entrez	16154
Expiration Date	Please enquire.